

Soluble Fms-Like Tyrosine Kinase 1 (sFlt-1) Antibody

| | |
|-------------------------|--|
| Catalog # | ABT-111009 |
| Clonality | Recombinant monoclonal |
| Host | Rabbit |
| Specificity | This antibody recognizes human sFlt-1. |
| Isotype | IgG |
| Cross Reactivity | NA |
| Concentration | 2.0 mg/ml (\pm 20%) |
| Buffer | PBS, pH 7.4 (8.1 mM Na ₂ HPO ₄ , 2.2 mM KH ₂ PO ₄ , 136.9 mM NaCl, 2.7 mM KCl) |
| Preservative | 0.05% Proclin 300 |
| Purity | \geq 90% |
| pI | NA |
| Applications | ELISA, Lateral Flow, Chemiluminescence Immunoassay |
| Storage | -20°C |
| Shelf Life | 36 months |

For research use only.