

N-terminal pro-B-type Natriuretic Peptide (NT-proBNP) Antibody

Catalog #	ABT-101044
Clonality	Recombinant monoclonal
Host	Sheep
Specificity	This antibody recognizes human NT-proBNP.
Isotype	IgG1
Cross Reactivity	NA
Concentration	2.0 mg/ml (\pm 20%)
Buffer	PBS, pH 7.4 (8.1 mM Na ₂ HPO ₄ , 2.2 mM KH ₂ PO ₄ , 136.9 mM NaCl, 2.7 mM KCl)
Preservative	0.05% Proclin 300
Purity	\geq 90%
pI	NA
Applications	ELISA, Lateral Flow, Chemiluminescence Immunoassay
Storage	-20°C
Shelf Life	36 months

For research use only.