

Recombinant human Renin protein (His tag)

Description

Product name	Recombinant human Renin protein (His tag)
Catalog#	ABT-9009
Known as	Angiotensinogenase
Expression system	HEK 293 cells
Tags	His tag C-terminus

Specifications

SDS-PAGE	50 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Lyophilized powder

Stability and Storage

Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Aliquot and store at -20°C for 12 months.
Constituents	Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 and 5% trehalose.


General information

Function	Renin is a highly specific endopeptidase, whose only known function is to generate angiotensin I from angiotensinogen in the plasma, initiating a cascade of reactions that produce an elevation of blood pressure and increased sodium retention by the kidney.
-----------------	--

Images

For research use only.

kDa M R



Human Renin (His Tag) on SDS-PAGE under reducing condition(R). The gel was visualized by Coomassie® Blue Staining. The purity of the protein is greater than 95%.

For research use only.