

Angiotensin I-BSA

Description

Product name Angiotensin I-BSA

Source Synthetic

Specifications

Format Liquid

Stability and Storage

Storage Store at -20°C, avoid repeated freeze/ thaw cycle.

Buffer 0.02M PBS, PH7.2

Stability. Aliquot and store at -20°C for 24 months.

General information

Function

Angiotensin is a kind of peptide substance with extremely strong effect of constricting blood vessel and stimulating adrenal cortex to secrete aldosterone. It participates in the regulation of blood pressure and body fluid. Angiotensin I is converted from renin action on angiotensinogen, and can reflect renin activity.

For research use only.