

## Testosterone-BSA

### Description

<b>Product name</b>	Testosterone-BSA
<b>Source</b>	Synthetic

### Specifications

<b>Format</b>	Liquid
---------------	--------

### Stability and Storage

<b>Storage</b>	Store at -20°C, avoid repeated freeze/ thaw cycle.
<b>Buffer</b>	0.02M PBS, PH7.2
<b>Stability.</b>	Aliquot and store at -20°C for 24 months.

### General information

<b>Function</b>	Testosterone is a steroid hormone secreted by the male testes or the female ovaries. It is the major male sex hormone, but it has important health effects for both men and women, including increased libido, strength, immune function, and fighting osteoporosis. Testosterone can be clinically used to diagnose primary testicular hypofunction, adrenal tumor, male precocious puberty and other diseases.
-----------------	--

**For research use only.**