

Recombinant FCys-C protein (His tag)

Description

Product name	Recombinant FCys-C protein (His tag)
Catalog#	ABT-9026
Known as	Cystatin C
Expression system	HEK 293 cells
Tags	His tag C-terminus

Specifications

SDS-PAGE	14.5 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Lyophilized powder


Stability and Storage

Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Aliquot and store at -20°C for 12 months.
Constituents	Lyophilized from 0.22 µm filtered solution in PBS, Ph 8.0.

General information

Function	Amyloid-beta binding; cysteine-type endopeptidase inhibitor activity; identical protein binding; protease binding.
-----------------	--

Images



FCys-C (His Tag) on SDS-PAGE under reducing condition(R). The gel was visualized by Coomassie® Blue Staining. The purity of the protein is greater than 95%.

For research use only.