

Recombinant Monkeypox Virus B14R protein (His tag)

Description

Product name	Recombinant Monkeypox Virus B14R protein (His tag)
Catalog#	ABT-9045
Known as	BR-209
Expression system	E.coli
Tags	His tag C-terminus

Specifications

SDS-PAGE	37 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Liquid


Stability and Storage

Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Aliquot and store at -20°C for 12 months.
Constituents	8M Urea.

General information

Function	IL-1 β -binding protein, with different fragments in Central African and West African strains, whose encoding gene has no crossover with variola virus.
-----------------	---

Images



Monkeypox Virus B14R (His Tag) on SDS-PAGE under reducing condition(R). The gel was visualized by Coomassie® Blue Staining. The purity of the protein is greater than 95%.

For research use only.