

Recombinant human SCCA protein (His tag)

Description

Product name	Recombinant human SCCA protein (His tag)
Catalog#	ABT-9108
Known as	Serpin B3; Protein T4-A; Squamous cell carcinoma antigen 1
Expression system	E.coli
Tags	His tag C-terminus

Specifications

SDS-PAGE	46 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Lyophilized powder

Stability and Storage


Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Store at -20°C for 12 months, or reconstitute for 3 months.
Constituents	Lyophilized from 0.22 µm filtered solution in Tris-HCl with 0.2M NaCl, pH8.

General information

Function	May act as a papain-like cysteine protease inhibitor to modulate the host immune response against tumor cells. Also functions as an inhibitor of UV-induced apoptosis via suppression of the activity of c-Jun NH2-terminal kinase (JNK1).
-----------------	--

Images

For research use only.



Human SCCA (His Tag) on SDS-PAGE under reducing condition(R). The gel was visualized by Coomassie® Blue Staining. The purity of the protein is greater than 95%.

For research use only.