

Recombinant FSAA protein (His tag)

Description

Product name	Recombinant FSAA protein (His tag)
Catalog#	ABT-9027
Known as	Serum amyloid A protein
Expression system	HEK 293 cells
Tags	His tag

Specifications

SDS-PAGE	13 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Lyophilized powder


Stability and Storage

Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Aliquot and store at -20°C for 12 months.
Constituents	Lyophilized from 0.22 µm filtered solution in PBS, Ph 7.4.

General information

Function	Major acute phase reactant. Apolipoprotein of the HDL complex.
-----------------	--

Images



FSAA (His tag) on SDS-PAGE under reducing condition(R).
The gel was visualized by Coomassie® Blue Staining. The
purity of the protein is greater than 95%.

For research use only.