

Recombinant Monkeypox Virus A44R protein (His tag)

Description

Product name	Recombinant Monkeypox Virus A44R protein (His tag)
Catalog#	ABT-9044
Known as	A44R
Expression system	E.coli
Tags	His tag C-terminus

Specifications

SDS-PAGE	15 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Liquid


Stability and Storage

Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Aliquot and store at -20°C for 12 months.
Constituents	PBS, pH7.4

General information

Function	Bind to the serum of sick monkeys, which has a significant IgM response, may be used as an antigen in the detection of monkeypox virus.
-----------------	---

Images



Monkeypox Virus A44R (His Tag) on SDS-PAGE under reducing condition(R). The gel was visualized by Coomassie® Blue Staining. The purity of the protein is greater than 95%.

For research use only.