

Recombinant human Hemopexin protein (His tag)

Description

Product name	Recombinant human Hemopexin protein (His tag)
Catalog#	ABT-9105
Known as	Beta-1B-glycoprotein; HPX
Expression system	HEK 293 cells
Tags	His tag C-terminus

Specifications

SDS-PAGE	37 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Lyophilized powder


Stability and Storage

Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Store at -20°C for 12 months, or reconstitute for 3 months.
Constituents	Lyophilized from 0.22 µm filtered solution in PBS, Ph7.4.

General information

Function	Binds heme and transports it to the liver for breakdown and iron recovery, after which the free hemopexin returns to the circulation.
-----------------	---

Images



Human Hemopexin (His Tag) on SDS-PAGE under reducing condition(R). The gel was visualized by Coomassie® Blue Staining. The purity of the protein is greater than 95%.

For research use only.