

Recombinant human AFP protein (His tag)

Description

Product name	Recombinant human AFP protein (His tag)
Catalog#	ABT-9005
Known as	Alpha-fetoprotein; Alpha-1-fetoprotein; Alpha-fetoglobulin
Expression system	HEK 293 cells
Tags	His tag C-terminus

Specifications

SDS-PAGE	70 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Lyophilized powder

Stability and Storage


Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Aliquot and store at -20°C for 12 months.
Constituents	0.22 µm filtered solution in PBS, pH7.4.

General information

Function	Binds copper, nickel, and fatty acids as well as, and bilirubin less well than, serum albumin. Only a small percentage (less than 2%) of the human AFP shows estrogen-binding properties.
-----------------	---

Images

For research use only.



Human AFP (His Tag) on SDS-PAGE under reducing condition(R). The gel was visualized by Coomassie® Blue Staining. The purity of the protein is greater than 95%.

For research use only.