

CG-BSA

Description

Product name	CG-BSA
Source	Synthetic

Specifications

Format	Liquid
---------------	--------

Stability and Storage

Storage	Store at -20°C, avoid repeated freeze/ thaw cycle.
Buffer	0.02M PBS, PH7.2
Stability.	Aliquot and store at -20°C for 24 months.

General information

Function	Cholyglycine (CG) is one of the conjugated cholic acids formed by the combination of cholic acid and glycine. CG is the most important bile acid component in late pregnancy serum. When hepatocytes are damaged, the ability of hepatocytes to ingest CG decreases, resulting in an increase in blood CG content.
-----------------	--

For research use only.