

E2-BSA

Description

Product name	E2-BSA
Source	Synthetic

Specifications

Format	Liquid
---------------	--------

Stability and Storage

Storage	Store at -20°C, avoid repeated freeze/ thaw cycle.
Buffer	0.02M PBS, PH7.2
Stability.	Aliquot and store at -20°C for 24 months.

General information

Function	<p>Estradiol (E2), also spelled oestradiol, is an estrogen steroid hormone and the major female sex hormone. It is involved in the regulation of the estrous and menstrual female reproductive cycles and is also one of the indicators to monitor follicular development and maturity. Examination of estradiol in blood and urine is of certain value for the diagnosis of endocrine and gynecological diseases such as sexual precocity and dysplasia.</p>
-----------------	---

For research use only.