

Recombinant human TNF- α protein (His tag)

Description

Product name	Recombinant human TNF- α protein (His tag)
Catalog#	ABT-9024
Known as	Prothrombin; Coagulation factor II
Expression system	HEK 293 cells
Tags	His tag C-terminus

Specifications

SDS-PAGE	17 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Lyophilized powder

Stability and Storage


Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Aliquot and store at -20°C for 12 months.
Constituents	Lyophilized from 0.22 μ m filtered solution in PBS, Ph 8.0.

General information

Function	Cytokine that binds to TNFRSF1A/TNFR1 and TNFRSF1B/TNFR. It is mainly secreted by macrophages and can induce cell death of certain tumor cell lines. Plays a role in angiogenesis by inducing VEGF production synergistically with IL1B and IL6.
-----------------	---

Images

For research use only.



Human TNF- α (His Tag) on SDS-PAGE under reducing condition(R). The gel was visualized by Coomassie® Blue Staining. The purity of the protein is greater than 95%.

For research use only.