

Recombinant human IL-5 protein (His tag)

Description

Product name	Recombinant human IL-5 protein (His tag)
Catalog#	ABT-9034
Known as	B-cell differentiation factor I; Eosinophil differentiation factor; T-cell replacing factor (TRF)
Expression system	HEK 293 cells
Tags	His tag C-terminus

Specifications

SDS-PAGE	16,18,20 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Lyophilized powder


Stability and Storage

Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Aliquot and store at -20°C for 12 months.
Constituents	Lyophilized from 0.22 µm filtered solution in PBS, pH7.4

General information

Function	Factor that induces terminal differentiation of late-developing B-cells to immunoglobulin secreting cells.
-----------------	--

Images



Human IL-5 (His Tag) on SDS-PAGE under reducing condition(R). The gel was visualized by Coomassie® Blue Staining. The purity of the protein is greater than 95%.

For research use only.



For research use only.

CONTACT US

**CHINA
GLOBAL**

400-100-1145
+1 312-802-6843

IVD@ptm-biolab.com
IVD@ptm-biolab.com

www.ptm-biolab.com.cn
www.ptmbio.com