

Recombinant human PCT protein (His tag)

Description

Product name	Recombinant human PCT protein (His tag)
Catalog#	ABT-9083
Known as	Calcitonin
Expression system	HEK 293 cells
Tags	His tag C-terminus

Specifications

SDS-PAGE	14.6 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Lyophilized powder

Stability and Storage

Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Store at -20°C for 12 months, or reconstitute for 3 months.
Constituents	Lyophilized from 0.22 µm filtered solution in PBS, pH7.4.


General information

Function	Calcitonin causes a rapid but short-lived drop in the level of calcium and phosphate in blood by promoting the incorporation of those ions in the bones; Katalcalcin is a potent plasma calcium-lowering peptide.
-----------------	---

Images

For research use only.

kDa M R



Human PCT protein on SDS-PAGE under reducing condition(R). The gel was visualized by Coomassie® Blue Staining. The purity of the protein is greater than 95%.

For research use only.