

Recombinant human P1NP protein (His tag)

Description

Product name	Recombinant human P1NP protein (His tag)
Catalog#	ABT-9030
Known as	Collagen alpha-1(I) chain; Alpha-1 type I collagen
Expression system	HEK 293 cells
Tags	His tag C-terminus

Specifications

SDS-PAGE	25 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Lyophilized powder


Stability and Storage

Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Aliquot and store at -20°C for 12 months.
Constituents	Lyophilized from 0.22 µm filtered solution in PBS, pH7.4

General information

Function	Type I collagen is a member of group I collagen (fibrillar forming collagen).
-----------------	---

Images



Human P1NP (His Tag) on SDS-PAGE under reducing condition(R). The gel was visualized by Coomassie® Blue Staining. The purity of the protein is greater than 95%.

For research use only.