

Recombinant human VEGF165 protein (His tag)

Description

Product name	Recombinant human VEGF165 protein (His tag)
Catalog#	ABT-9039
Known as	Vascular endothelial growth factor A
Expression system	HEK 293 cells
Tags	His tag C-terminus

Specifications

SDS-PAGE	21 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Lyophilized powder

Stability and Storage


Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Aliquot and store at -20°C for 12 months.
Constituents	Lyophilized from 0.22 µm filtered solution in PBS, pH7.4

General information

Function	Growth factor active in angiogenesis, vasculogenesis and endothelial cell growth. Induces endothelial cell proliferation, promotes cell migration, inhibits apoptosis and induces permeabilization of blood vessels.
-----------------	--

Images

For research use only.



Human human VEGF 165 (His Tag) on SDS-PAGE under reducing condition(R). The gel was visualized Coomassie® Blue Staining. The purity of the protein is greater than 95%.

For research use only.