

Recombinant human TSH protein (His tag)

Description

Product name	Recombinant human TSH protein (His tag)
Catalog#	ABT-9071
Known as	<i>Thyroid-stimulating hormone</i>
Expression system	HEK 293 cells
Tags	His tag C-terminus

Specifications

SDS-PAGE	26.3 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Lyophilized powder

Stability and Storage


Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Store at -20°C for 12 months, or reconstitute for 3 months.
Constituents	Lyophilized from 0.22 µm filtered solution in PBS, pH7.4.

General information

Function	Thyrotropin (TSH) is a thyroid-stimulating hormone produced by thyrotrope cells in the anterior pituitary gland. Thyrotropin regulates the endocrine function of the thyroid.
-----------------	---

Images

For research use only.



Human TSH (His Tag) on SDS-PAGE under reducing condition(R). The gel was visualized by Coomassie® Blue Staining. The purity of the protein is greater than 95%.

For research use only.