

Recombinant FeLV protein (His tag)

Description

Product name	Recombinant FeLV protein (His tag)
Catalog#	ABT-9043
Known as	Gag polyprotein; Core polyprotein
Expression system	E.coli
Tags	His tag C-terminus

Specifications

SDS-PAGE	26 kDa, reducing conditions
Purity	>95% SDS-PAGE
Form	Lyophilized powder

Stability and Storage


Storage	Store at -20 ~ -80°C, avoid repeated freeze/ thaw cycle
Stability	Aliquot and store at -20°C for 12 months.
Constituents	Lyophilized from 0.22 µm filtered solution in PBS, Ph 8.0.

General information

Function	Plays a role in budding and is processed by the viral protease during virion maturation outside the cell. During budding, it recruits, in a PPXY-dependent or independent manner, Nedd4-like ubiquitin ligases that conjugate ubiquitin molecules to Gag, or to Gag binding host factors.
-----------------	---

Images

For research use only.



FeLV (His Tag) on SDS-PAGE under reducing condition(R). The gel was visualized by Coomassie® Blue Staining. The purity of the protein is greater than 95%.

For research use only.