

## Progesterone-BSA

### Description

<b>Product name</b>	Progesterone-BSA
<b>Source</b>	Synthetic

### Specifications

<b>Format</b>	Liquid
---------------	--------

### Stability and Storage

<b>Storage</b>	Store at -20°C, avoid repeated freeze/ thaw cycle.
<b>Buffer</b>	0.02M PBS, PH7.2
<b>Stability.</b>	Aliquot and store at -20°C for 24 months.

### General information

#### Function

Progesterone secreted by the corpus luteum of the ovary, has a significant morphological effect on endometrium and is necessary for the maintenance of pregnancy. Progesterone content was closely related to ovarian corpus luteum and pregnant placenta. It is clinically used for determining ovulation, detecting corpus luteum function and pregnancy progression, and for the diagnosis of obstetrical diseases such as threatened abortion, habitual abortion and ectopic pregnancy.

**For research use only.**