

Total TAU Protein (T-TAU) Antibody

Catalog #	ABT-112008
Clonality	Recombinant monoclonal
Host	Rabbit
Specificity	This antibody recognizes human T-TAU.
Isotype	IgG
Cross Reactivity	NA
Concentration	2.0 mg/ml (\pm 20%)
Buffer	37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl
Preservative	0.05% Proclin 300
Purity	\geq 90%
pI	NA
Applications	ELISA, Lateral Flow, Chemiluminescence Immunoassay
Storage	-20°C
Shelf Life	36 months

For research use only.

CONTACT US

CHINA

400-100-1145

IVD@ptm-biolab.com

www.ptm-biolab.com.cn

GLOBAL

+1 312-802-6843

info@ptmbio.com

www.ptmbio.com