

p-Tau 217 Antibody

Catalog #	ABT-112013
Clonality	Recombinant monoclonal
Host	Rabbit
Specificity	This antibody recognizes human p-Tau 217.
Isotype	IgG
Cross Reactivity	NA
Concentration	2.0 mg/ml (\pm 20%)
Buffer	37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl
Preservative	0.05% Proclin 300
Purity	\geq 90%
pI	NA
Applications	ELISA, Lateral Flow, Chemiluminescence Immunoassay
Storage	-20°C
Shelf Life	36 months

For research use only.