

Total TAU Protein (T-TAU) Antibody

| | |
|-------------------------|---|
| Catalog # | ABT-112009 |
| Clonality | Recombinant monoclonal |
| Host | Mouse |
| Specificity | This antibody recognizes human T-TAU. |
| Isotype | IgG1 |
| Cross Reactivity | NA |
| Concentration | 2.0 mg/ml (\pm 20%) |
| Buffer | 37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl |
| Preservative | 0.05% Proclin 300 |
| Purity | \geq 90% |
| pI | NA |
| Applications | ELISA, Lateral Flow, Chemiluminescence Immunoassay |
| Storage | -20°C |
| Shelf Life | 36 months |

For research use only.

CONTACT US

CHINA

400-100-1145

IVD@ptm-biolab.com

www.ptm-biolab.com.cn

GLOBAL

+1 312-802-6843

info@ptmbio.com

www.ptmbio.com