



Product Datasheet

Name	Neuron Specific Enolase (NSE) Antibody
Catalog#	ABT-102071
Format	Monoclonal antibody, in vitro recombination
Host / Source	Mouse
Specificity	Identify Human NSE
Class / isotype	IgG1
Cross-reactivity	NA
Concentration	2.0 mg/ml (\pm 20%)
Supplied in	PBS pH 7.4 Buffer Solution (8.1 mM Na ₂ HPO ₄ , 2.2 mM KH ₂ PO ₄ , 136.9 mM NaCl, 2.7 mM KCl)
Preservative	0.05% Proclin 300
Appearance	Transparent clear. Precipitation may occur after long-term storage
Purity	\geq 90%
pI	NA
Applications	Elisa, Lateral flow, Chemiluminescence

Matched pairs

		Detection					
		102066	102067	102068	102069	102070	102071
Capture	102066	-	+	+	+	+	+
	102067	+	-	+	+	+	+
	102068	+	+	-	+	+	+
	102069	-	+	+	-	+	+
	102070	+	+	+	+	-	+
	102071	+	+	+	+	+	-

Storage -20 °C

Expiration 36 months

For research use only.